

Your perfect sized condom!

MY.SIZE

Find your perfect sized My.Size condom for a more naturally safe sensation.

My.Size is the customised condom brand that come in seven sizes ranging from 47mm to 69mm. With My.Size, you will always find the perfect fitting condom for a more naturally safe sensation.



Cut out Measuring Tape

Cut along the dotted line of the measuring tape with a pair of scissors.


Be aware that you can get sharp edges from cutting paper.



Measure

The penis circumference is determined by taking the measuring tape and wrapping it around the penis. The result you can read on the lines marked on the tape.

The **penis circumference** is vital for choosing the right size of My.Size condom.




Establishing your My.Size size

In the chart, you will find your correct My.Size size.

For example: If the circumference of your penis is 11.5cm, then your ideal My.Size size is 53.

If you wish to have a tighter sitting condom or you would like to have a greater sensation, choose a smaller or bigger size from the My.Size condom range.

! IMPORTANT: When printing this page, click „None“ for the Page Scaling option. Otherwise the measuring tape will not be the correct size.



MY.SIZE size **Penis circumference** in cm

47	9,5 – 10,0
49	10,0 – 11,0
53	11,0 – 11,5
57	11,5 – 12,0
60	12,0 – 13,0
64	13,0 – 14,0
69	14,0 – 15,0



For more information visit www.measureyourpenis.today